

	Patient	Controls
	Counts	Count range
Neutrophils (x10 <sup>3</sup> /ml)	4200	1500-7400
Eosinophils (x10 <sup>6</sup> /ml)	100	50-1290
Basophils (x10 <sup>3</sup> /ml)	100	0-240
Monocytes (x10 <sup>3</sup> /ml)	700	330-1240
	Counts	Counts: mean (range)
cDC1 (/ml)	438	860 (0-2160)
cDC2 (/ml)	14800	15600 (6580-32800)
pDC (/ml)	23800	11200 (4510-20400)
	Frequency	Frequency: mean (range)
cDC1 (%)	0.0044	0.02 (0-0.04)
cDC2 (%)	0.148	0.27 (0.09-0.42)
pDC (%)	0.238	0.19 (0.09-0.37)

Supplemental Table 3: Myeloid cell counts. Counts are expressed per ml of blood. cDC1: conventional dendritic cells type 1 (CD11c<sup>+</sup>CD141<sup>+</sup>); cDC2: conventional dendritic cells type 2 (CD11c<sup>+</sup>CD1c<sup>+</sup>); pDC: plasmacytoid dendritic cells (CD303<sup>+</sup>).